
Disposable Needle Guide

For use with Toshiba PVM-740RT
and PVT-770RT transducers

Reference Guide

Precautions for use

English

WARNING

- When PVM-740RT / PVT-770RT transducer is connected to ultrasound system, guideline for Toshiba biopsy adapter UAGV-023A is displayed. Needle angle of CIVCO needle guidance system differs from UAGV-023A. To use CIVCO needle guidance system, optional software must be installed. After installation, guideline for UAGV-023A is not displayed. Contact your Toshiba representative for details.
- The dedicated software must be installed to display guideline for CIVCO needle guidance system. Using CIVCO needle guidance system with UAGV-023A guideline, biopsy may be performed at an inaccurate biopsy angle, resulting in serious injury to patient.
- On-screen guidelines are a reference for needle direction. Guidelines do not indicate actual position of needle.

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

INTENDED USE

- This device provides physicians with a tool for performing needle/instrument guided procedures with the use of the diagnostic ultrasound endocavity transducer.

WARNING

- *In the USA, federal law restricts device to sale or use by or on order of a physician.*
- *Before use, you should be trained in ultrasonography. For instructions on the use of your transducer, see your system's manual.*
- *Transducer covers containing natural rubber latex may cause allergic reactions. In the USA, refer to FDA Medical Alert MDA91-1.*
- *Disposable components are packaged sterile and are single-use only. Do not use if integrity of packaging is violated or if expiration date has passed.*
- *For illustration purposes only, transducer may be shown without a transducer cover. Always place a cover over transducer to protect patients and users from cross-contamination.*
- *Use only water-soluble agents or gels. Petroleum or mineral oil based materials may harm cover.*
- *Do not reuse, reprocess or resterilize single-use device. Reuse, reprocessing or resterilization may create a risk of contamination of the device, cause patient infection or cross-infection.*

I. Attaching Bracket To Transducer

WARNING

- *Bracket is packaged non-sterile and is reusable. To avoid possible patient contamination, ensure bracket is properly cleaned, sterilized or disinfected before each use.*
- *If bracket does not fit properly, do not use.*

1. Attach bracket to transducer by aligning locating features. Ensure bracket is firmly attached.

II. Covering The Transducer

1. Place an appropriate amount of gel inside cover and/or on transducer face. Poor imaging may result if no gel is used.
2. Insert transducer into cover making sure to use proper sterile technique. Pull cover tightly over transducer face to remove wrinkles and air bubbles, taking care to avoid puncturing cover.
3. Secure cover with enclosed bands.
4. Inspect cover to ensure there are no holes or tears.

III. Attaching Needle Guide To Bracket

1. Using proper sterile technique, align guide with locating features of transducer and bracket.
2. Place front of needle guide onto front of transducer and snap back end into slot on bracket.

IV. Needle Path Verification

WARNING

- *Before using component for first time, perform needle path verification to verify system and biopsy guide relationships.*
- *Use appropriate needle length to reach target area.*
- *See transducer manual for proper procedure.*
- *Use care when inserting needle to avoid bending.*

V. Disposal

WARNING

- *Dispose of single-use components as infectious waste.*
- *Users of this product have an obligation and responsibility to provide the highest degree of infection control to patients, co-workers and themselves. To avoid cross-contamination, follow infection control policies established by your facility.*
- *See your system's manual for reprocessing transducer between use.*

WARNING**VI. Cleaning, Disinfecting and Sterilizing Bracket**

- Users of this product have an obligation and responsibility to provide the highest degree of infection control to patients, co-workers and themselves. To avoid cross-contamination, follow infection control policies established by your facility.
- See your system's user guide for reprocessing transducer between use.
- The recommended protocol for reprocessing the reusable bracket has been evaluated for effectiveness and compatibility. If choosing to use a reprocessing protocol other than those listed, a validation must first be completed to ensure effectiveness and compatibility.
- Do not gas sterilize or autoclave bracket.

CLEANING THE BRACKET

1. After each use, remove bracket and needle guide from transducer, discarding single-use, disposable components. Remove visible contaminants from bracket surface using a small, soft instrument brush. Keep bracket from drying out until complete cleaning can be accomplished.
2. Soak bracket for minimum of five minutes in neutral pH, low foaming enzymatic detergent, such as ENZOL[®] Enzymatic Detergent (Johnson & Johnson). While immersed, use instrument brush to remove trapped contaminants from surfaces and features. If visible contaminants cannot be easily removed, repeat soaking procedure for an additional five minutes. Remove bracket from cleaning solution and remove any remaining residue with dry wipe. Follow cleaning solution manufacturer's directions for use and recommendations for concentration.

DISINFECTING AND STERILIZING BRACKET

- High-level disinfect bracket using CIDEX OPA[®] ortho-Phthalaldehyde Solution (Johnson & Johnson) or equivalent .55% ortho-phthalaldehyde-based solution. Follow manufacturer's instructions and recommendations for concentration, time of contact and post-process procedure.
- High-level disinfect or sterilize bracket using CIDEX[®] Activated Dialdehyde Solution (Johnson & Johnson) or equivalent 2% glutaraldehyde-based solution, CIDEX Plus[®] (Johnson & Johnson) or equivalent 3.4% glutaraldehyde-based solution, or a hydrogen peroxide-based solution. Follow manufacturer's instructions and recommendations for concentration, time of contact and post-process procedure.
- Inspect bracket for damage such as cracks or breakage. If damage is evident, discontinue use of bracket and contact your ultrasound system support representative.

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COPYRIGHT © 2012 ALL RIGHTS RESERVED. CIVCO IS A REGISTERED TRADEMARK OF CIVCO MEDICAL SOLUTIONS. ALL OTHER TRADEMARKS ARE PROPERTY OF THEIR RESPECTIVE OWNERS. PRINTED IN USA. IN USA: ALL PRODUCTS MAY NOT BE LICENSED IN ACCORDANCE WITH CANADIAN LAW.
MPS GmbH